## **Trimethoprim**

# 80 mg in Co-Formulation Tablets with 400 mg Sulfamethoxazole

(A screening method for sulfamethoxazole in this combined formulation is available in GPHF Minilab Manual Volume II, Method 6.39, pp. 184-187)

#### **Structure:**

**Molecular Formula and Mass:** C<sub>14</sub>H<sub>18</sub>N<sub>4</sub>O<sub>3</sub> - 290.323

Category: Antibacterial

### Sample:

Grind one tablet and dissolve in 100 mL of methanol. Shake for at least 10 min and filter. Dilute 1 mL of the stock solution with 2 mL of methanol. Final concentration of sample solutions = 0.267 mg/mL, which is the required concentration representing 100%.

## **Standards:**

#### High Standard:

The high limit is 115%; therefore the concentration of the high standard =  $0.267 \text{ mg/mL} \times 1.15 = 0.307 \text{ mg/mL}$ . Weigh approximately 30.7 mg of standard and dissolve it in 100 mL methanol. If you weighed 30.6 mg of standard, dissolve it in:  $30.6 \text{ mg} \div 0.307 \text{ mg/mL} = 99.7 \text{ mL}$  of methanol. This makes the high standard solution concentration equal to 0.307 mg/mL. Low Standard:

The low limit is 85%; therefore the concentration of the low standard =  $(0.267 \text{ mg/mL} \times 0.85 = 0.227 \text{ mg/mL}$ . Dilute 3.00 mL of high standard to 4.00 mL by adding 1.00 mL of methanol. This gives a concentration of 0.307 mg/mL  $\times$  3.00 mL  $\div$  4.00 mL = 0.230 mg/mL, which is 86.2%.

## **Spotting:**

Spot on the 5  $\times$  10 cm silica gel TLC aluminum plate with 3.00  $\mu$ L aliquots as follows:

Left spot low standard (85%) =  $0.690 \mu g$ 

Center Spot 100% sample =  $0.801 \mu g$ 

Right Spot high standard (115%) =  $0.921 \mu g$ 


#### **Development:**

Mix 24.0 mL of ethyl acetate and 8.00 mL of methanol. Develop the plate in a small glass chamber with approximately 20.0 mL of this solution until the solvent front reaches within 1 cm of the top of the TLC plate.

 $(R_f = 0.21)$ 

#### **Detection:**

<u>UV:</u> Dry the plate and observe under ultraviolet light at 254 nm. Observe the intensities and the sizes of the spots.



Developed and tested by Bingsong Zeng and Joseph Sherma
Department of Chemistry, Lafayette College, Easton, PA, USA
November 2017

G's research was supported by the Lafayette College EXCEL Scholars I

Bingsong Zeng's research was supported by the Lafayette College EXCEL Scholars Program